

Disposable Biopsy Valves

Alkapharm Disposable Biopsy Valves are used to cover the openings to the biopsy/suction channel inlets of Olympus®, Fujinon® or Pentax® gastrointestinal endoscopes.

The biopsy valves provide access for endoscopic device passage and exchange, help to maintain insufflation, and minimize the leakage of bioburden from the biopsy port throughout the gastrointestinal endoscopic procedure. The disposable biopsy valves are single use which eliminates the need for manual cleaning and reprocessing. Disposable biopsy valves reduce cross contamination helping to prevent potential patient/staff safety risks.

FEATURES

- No “splitting” or leakage for general GI procedures
- Helps reduce cross-contamination
- Eliminates the need for reprocessing dirty valves
- Latex free
- Disposable, single use
- Individually bagged, available in boxes of 100 units.
- CE Marked
- Distributed by Mediteq Healthcare Solutions, available through United Drug.
 - Product code (ALK-BVO) - Olympus & Fujinon
 - Product code (ALK-BVP) - Pentax

