

Neutra-Cleanse



Non foaming neutral pH detergent for cleaning of surgical instruments & endoscopes

- Super-concentrated formula, 2-3mls / litre
- 5L container will make up 2,500 litres of ready-to-use solution.
- Solvent Free
- Added water softener
- Neutral pH: pH diluted in hard water 7.0 +/- 0.5
- Fragrance free
- Colour: Blue
- No residues & free rinsing
- Enzyme free
- Compatible & safe to use with all makes of flexible endoscopes
- **Why neutral?**
 - A review from the Health & Safety Authority in England drove the change from Enzymatic detergents to Neutral detergents for instrument & endoscope pre-cleaning in England, Scotland & Wales.
 - The BSG 2014 guidelines also acknowledges this move away from enzymatics, see below;
 - *“Previous editions of this guidance have recommended the use of enzymatic detergents based on their theoretical ability to digest mucus and other biological material from within narrow endoscope channel lumens. In the CFPP 01- 06 guidance for health service units in England no preference for enzymatic detergents is expressed (19). A move away from enzymatic detergents has been prompted by reports of occupational asthma and skin sensitization (20).”*

New !

