

# Disposable Isolation Gowns



Disposable isolation gowns are used to protect healthcare workers exposed body areas and prevent contamination of clothing with blood, body substances and other potentially infectious material.\*

Meditaq's disposable isolation gowns are impervious / fluid resistant, but also breathable, designed with healthcare workers safety and comfort in mind.

These disposable isolation gowns are CE marked and manufactured in an ISO Accredited facility. They also come with wall mounted dispensers for ease of access.



## SPECIFICATIONS

<b>Classification</b>	Personal Protective Equipment (PPE) cat. III
<b>Certification</b>	<ul style="list-style-type: none"> <li>- Manufactured according to ISO 13485 standards</li> <li>- CE-certified (notified body 0624)</li> <li>- EC Declaration of Conformity according to EU Regulation 2016/425</li> <li>- PPE cat. III certified (EU Type Examination certificate)</li> </ul>
<b>Material</b>	Polypropylene (PP) non-woven & Polyethylene (PE)
<b>Thickness</b>	35gsm: 23gsm PP + 12gsm PE
<b>Dimensions</b>	120cm x 140cm
<b>Cuffs</b>	Knitted Polyester
<b>Norms</b>	<ul style="list-style-type: none"> <li>- EN ISO 13688:2013+A1:2021</li> <li>- EN 14126:2003+AC:2004</li> <li>- EN 14325:2004</li> </ul>
<b>Packaging information</b>	Packed per 10 pieces, 10 packs of 10 pieces per carton
<b>Additional specifications</b>	Latex Free, Non-sterile & Single-use only

Code	Product	Description
MTQ-GL-Yellow	35gsm IsolationGown Yellow	100 pack, PE+ PP Non-woven/Knitted Cuff/Universal Size

For more information, please call **01 816 2015**  
Email: [info@mediteq.ie](mailto:info@mediteq.ie) or visit [www.mediteq.ie](http://www.mediteq.ie)



All of our products are available to order through United Drug Distribution.  
For more information please contact your Mediteq account manager.

